

Our BESSEY tools in use

Repair of a cabinet door



Repair of doors from a solid wood cabinet

- The mitred and glued door frames had come away from the glue in places and the centre panel had shifted completely out of place. The frame must be re-glued.
- Two K Body REVO clamps were set, aligned on the bench. The frame was then placed in and two further K Body REVO clamps placed over the top.
- Then all of the mitres could be very simply and precisely pressed together by careful turning of the four clamps, without the frames being pressed out of true in the process.

Which BESSEY products were used and what characterises them:

Four **K Body REVO KR** clamps with an opening of 800 mm were used.

- These are perfectly suited to sensitive gluing and clamping of mitres.
- There is generally no need for additional supports here and this makes the re-gluing much simpler.
- Gluing right-angled joints is effortless with the body clamps.



These BESSEY products were used:



K Body REVO KR

- Clamping force up to 7,000 N
- Very large parallel clamping surface
- With 3 pressure caps to protect delicate surfaces
- Can be converted to spreading function with no need for tools
- High-quality 2-component plastic handle

Further products from the K Body REVO range:



Vario K Body REVO KRV



Framing set KP



Tilting K Body clamp adapter KR-AS



K Body Extender KBX

Contact information:

Master joiner / author
Heiko Rech
Spiemontstrasse 38
66 606 St. Wendel
Germany

Telephone: +49 163 - 60 41 947
Email: info@heiko-rech.de
Website: <http://holzwerkerblog.de>

For information about
BESSEY products, please visit
www.bessey.de